Flo-Rester®

Disposable Internal Vessel Occluder

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SYMBOL DEFINITIONS:

The contents of the unopened, undamaged pouches are sterile. Flo-Rester is Ethylene Oxide sterilized.

This product is for single use only, do not resterilize.

The date by which the product must be used.

Lot number

Catalog number

Read instructions for use prior to using this product.

Shaft length

Bulb diameter

PN: SSI part number. TN: SSI tracking number. IC: SSI internal code.

INDICATIONS:

Flo-Resters® are used to internally occlude blood vessels during anastomosis. The occluder bulbs are inserted proximally and distally through the vessel opening. The occluder bulbs allow a dry field to be maintained during the anastomosis of a bypassed graft or the suturing of two cut ends of a vessel. A pull tab allows the Flo-Rester to be easily removed prior to the final sutures being tied. Select the Flo-Rester size which is appropriate for the internal diameter of the blood vessel.

DESCRIPTION:


Flo-Rester is a sterile, single use, disposable occluder composed of two silicone rubber bulbs on the distal ends of a flexible shaft. A tab on a tether is attached to the midpoint of the flexible shaft forming a “T” configuration. Flo-Rester bulbs and tab are radiopaque.

INSTRUCTIONS FOR USE:

Each surgeon is responsible for the proper procedure and technique to be used with this device.

The Flo-Rester size is chosen based upon the internal diameter of the vessel.

The flexible shaft is held up by pickups adjacent to the bulb being inserted. The Flo-Rester is designed to easily slide inside the vessel. Never force the Flo-Rester into the vessel as subintimal dissection could occur. The first bulb is inserted through the vessel opening and into one side of the vessel. The second bulb is inserted in a similar manner. When correctly placed, the Flo-Rester flexible shaft will be straight. The pull tab is laid to one side out of the way (see illustration below). Unless there are perforating vessels between the bulbs of the Flo-Rester, bleeding from the vessel is completely controlled.

The Flo-Rester allows the blood vessels to be easily visualized since the walls of the vessels are separated by the bulbs.

WARNINGS:

- Due to the possible danger of tearing because of the friability of some vessels, care must be used in dealing with these vessels.
- Sub-intimal dissection must be guarded against. NEVER force the Flo-Rester bulb into a vessel.
- Carefully choose the correct Flo-Rester size according to the bulb diameter required.
- During the anastomosis of the graft to the artery, care must be exercised so the suture needle will not pierce the Flo-Rester's bulb or flexible shaft. The silicone rubber is soft and similar to tissue in consistency. Thus, it is not easy to feel the needle entering the bulb or shaft.

- Do not resterilize. The Flo-Rester is a disposable, single use device and is not intended for reuse.

DISCLAIMER OF WARRANTIES:

Synovis Surgical Innovations (SSI), a division of Synovis Life Technologies, Inc., warrants that reasonable care has been used in the manufacture of this device. This warranty is exclusive and in lieu of all other warranties whether expressed, implied, written or oral, including, but not limited to, any implied warranties of merchantability or fitness. As a result of biological differences in individuals, no product is 100% effective under all circumstances. Because of this fact and since SSI has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after it leaves its possession, SSI does not warrant either a good effect or against an ill effect following its use. The manufacturer shall not be liable for any incidental or consequential loss, damage or expense arising directly or indirectly from the use of this device. SSI will replace any device which is defective at the time of shipment. No representative of SSI may change any of the foregoing or assume any additional liability or responsibility in connection with this device.

Rx ONLY CAUTION:

Federal (USA) law restricts this device to sale by, or on the order of, a physician.

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Flo-Rester® is a registered trademark of Synovis Life Technologies, Inc.
Flo-Rester® ist ein eingetragenes Warenzeichen der Synovis Life Technologies, Inc.
Flo-Rester® est une marque de fabrique de Synovis Life Technologies, Inc.
Flo-Rester® è un marchio depositato della Synovis Life Technologies, Inc.
Flo-Rester® es una marca registrada de Synovis Life Technologies, Inc.
Flo-Rester® is een gedeponeerd handelsmerk van Synovis Life Technologies, Inc.