Synovis™
Surgical Innovations

Flo-Rester®

Reusable Internal Vessel Occluder
Interner Gefäßokkluder für den mehrfachen Gebrauch
Obturateur vasculaire interne réutilisable
Dispositivo riutilizzabile per l'occlusione interna del vaso
Ocluser vascular interno reutilizable
Herbruikbare inwendige bloedvatsluiser
Til intern karokkludering
Återanvändbar intern kärlockkluderare

INSTRUCTIONS FOR USE	2
GEBRAUCHSANWEISUNG	5
INSTRUCTIONS D'UTILISATION	8
ISTRUZIONI PER L'USO	11
INSTRUCCIONES DE USO	14
GEBRUIKSAANWIJZING	17
BRUGSANVISNING	20
BRUKSANVISNING	23

CE
0086

SYMBOL DEFINITIONS:

Can not be used more than 5 times

The date by which the product must be used

Lot number

Catalog number

Read instructions for use prior to using this product.

Shaft length

Bulb diameter

PN:SSI part number. TN:SSI tracking number. IC:SSI internal code.

INDICATIONS

Flo-Resters are used to internally occlude blood vessels during anastomosis. The occluder bulbs are inserted proximally and distally through the vessel opening. The occluder bulbs allow a dry field to be maintained during the anastomosis of a bypass graft or the suturing of two cut ends of a vessel. A pull tab allows the Flo-Rester to be easily removed prior to the final sutures being tied. Select the Flo-Rester size which is appropriate for the internal diameter of the blood vessel.

DESCRIPTION

Flo-Resters are composed of two silicone rubber bulbs on the distal ends of a flexible shaft. A tab on a tether is attached to the midpoint of the flexible

shaft, thus forming a "T" configuration. The Flo-Rester bulbs and pull tabs are radiopaque.

INSTRUCTIONS FOR USE

Each surgeon is responsible for the proper procedure and technique to be used with this device.

Before using the Flo-Rester, the bulbs must be carefully inspected to insure secure attachment to the flexible shaft and that the Flo-Rester is in satisfactory condition. Any questionable device must not be used.


I. HANDLING

Silicone elastomers used in the production of this product are among the most nonreactive materials known for this application. Great care has been taken to keep them free from contaminants such as dust, lint, oil, etc. However, silicone attracts these contaminants readily; consequently the package of the enclosed product should be opened in as clean an atmosphere as can be provided.

II. USAGE

The Flo-Rester size is chosen based upon the internal diameter of the vessel. The flexible shaft is held by pickups adjacent to the bulb to be inserted. The Flo-Rester is designed to easily slide inside the vessel. Never force the Flo-Rester into the vessel as subintimal dissection could occur. The first bulb is inserted through the vessel opening and into one side of

the vessel. The second bulb is inserted in a similar manner. When correctly placed the Flo-Rester flexible shaft will be straight. The pull tab is laid to one side out of the way (see illustration below). Unless there are perforating vessels between the bulbs of the Flo-Rester, bleeding from the vessel is completely controlled.

The Flo-Rester allows the blood vessels to be easily visualized since the walls of the vessels are separated by the bulbs.

III. CLEANING

After use, the Flo-Rester should be carefully washed with enzymatic detergent according to detergent manufacturer's instructions. All residue should be removed. A final rinse in non-pyrogenic water is recommended.

If blood is seen in the tubing, it indicates that the tubing has been punctured by a needle or similar object. When this occurs, the Flo-Rester should be replaced as it is impossible to remove the blood from the tubing.

IV. STERILIZATION

The Flo-Rester must be removed from the pouch and then should be wrapped in a clean, lint-free surgical towel or other suitable lint-free wrapping. It is recommended that the device be sterilized **NO MORE THAN FIVE (5) TIMES** by one of the following methods:

- Ethylene Oxide (EtO) gas sterilization according to standard validated hospital procedures
- Steam sterilization at 134°C for 3.5 minutes (270°F, U.S.)
- Steam sterilization at 121°C for 15 minutes (250°F, U.S.)

WARNINGS

- Flo-Resters must be sterilized before use (see STERILIZATION).
- Due to the possible danger of tearing because of the friability of some vessels, care must be used in dealing with these vessels.
- Sub-intimal dissection must be guarded against. NEVER force the Flo-Rester bulb into a vessel.
- Carefully choose the correct Flo-Rester size according to bulb diameter required.
- During anastomosis of the graft to the artery, care must be exercised so the suture needle will not pierce the Flo-Rester's bulb or shaft. The silicone rubber is soft and similar to tissue in consistency. Thus, it is not easy to feel the needle entering the bulb or shaft.

- Care must be exercised in cleaning the Flo-Rester (see CLEANING). Flo-Rester must be inspected prior to use to insure the bulbs remain securely fastened to the shaft and to ensure the device is in satisfactory condition.

PRECAUTIONS

Before using the Flo-Rester, the bulbs must be carefully inspected to insure secure attachment to the flexible shaft and that the Flo-Rester is in satisfactory condition. Any questionable device must not be used.

DISCLAIMER OF WARRANTIES:

Synovis Surgical Innovations (SSI), a division of Synovis Life Technologies, Inc., warrants that reasonable care has been used in the manufacture of this device. This warranty is exclusive and in lieu of all other warranties whether expressed, implied, written or oral, including, but not limited to, any implied warranties of merchantability or fitness. As a result of biological differences in individuals, no product is 100% effective under all circumstances. Because of this fact and since SSI has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after it leaves its possession, SSI does not warrant either a good effect or against an ill effect following its use. The manufacturer shall not be liable for any incidental or consequential loss, damage or expense arising directly or indirectly from the use of this device. SSI will replace any device which is

defective at the time of shipment. No representative of SSI may change any of the foregoing or assume any additional liability or responsibility in connection with this device.

Rx ONLY CAUTION:

Federal (USA) law restricts this device to sale by, or on the order of, a physician.

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Flo-Rester® is a registered trademark of Synovis Life Technologies, Inc.
Flo-Rester® ist ein eingetragenes Warenzeichen der Synovis Life Technologies, Inc.
Flo-Rester® est une marque de fabrique de Synovis Life Technologies, Inc.
Flo-Rester® è un marchio depositato della Synovis Life Technologies, Inc.
Flo-Rester® es una marca registrada de Synovis Life Technologies, Inc.
Flo-Rester® is een gedeponeerd handelsmerk van Synovis Life Technologies, Inc.