Synovis™
Surgical Innovations

Flo-Thru Intraluminal Shunt®

Intraluminal Shunt
Intraluminal-Shunt
Dérivation intraluminale
Shunt intraluminale

Derivación intraluminal
Intraluminale shunt
Intraluminal shunt
Intraluminal shunt

INSTRUCTIONS FOR USE	2
GEBRAUCHSANWEISUNG	4
INSTRUCTIONS D'UTILISATION	7
ISTRUZIONI PER L'USO	10
INSTRUCCIONES DE USO	13
GEBRUIKSAANWIJZING	16
BRUGSANVISNING	19
BRUKSANVISNING	22

CE
0086

SYMBOL DEFINITIONS

	The contents of the unopened, undamaged pouches are sterile. Flo-Thru is sterilized with Ethylene Oxide.
	This product is for single use only, do not resterilize.
	The date by which the product must be used.
	Lot number
	Catalog number
	Read instructions for use prior to using this product.
	Overall length
	Bulb diameter

PN:SSI part number. TN:SSI tracking number. IC:SSI internal code.

INDICATIONS

The Flo-Thru Intraluminal Shunt is indicated for use in coronary artery or peripheral vascular procedures. The device shunts blood at the anastomotic site which provides a temporary blood free operative field for suturing while allowing blood to flow distal to the anastomosis. The device is removed just prior to the final suturing of the vessel.

DESCRIPTION


The sterile, single use Flo-Thru Intraluminal Shunt is a one-piece radiopaque silicone tube with atraumatic bulbs shaped at each end.

Openings at the ends of the bulbs allow blood to flow through the shunt and distal to the anastomosis. A radiopaque tab, which identifies the outer diameter of the bulbs, is attached to the shunt by a tether to facilitate positioning and removal of the shunt.

INSTRUCTIONS FOR USE

Each surgeon is responsible for the proper procedure and technique to be used with this device.

- Select the Flo-Thru Intraluminal Shunt according to the internal diameter of the vessel.
- Hold the shunt with a forceps adjacent to the bulb to be inserted.
- Insert the first bulb into the proximal side of the arteriotomy. A flush of blood from the opposite end of the shunt confirms the intraluminal position and facilitates the removal of air. The other bulb can then be inserted into the arteriotomy.
- The shunt is designed to easily slide inside the vessel. Never force the shunt into the vessel as intimal dissection could occur.
- When correctly placed, the shunt will be straight. Place the tab to one side, out of the way (see illustration).
- Unless there are perforating vessels between the bulbs of the shunt and if the appropriate size has been selected, bleeding from the vessel is completely controlled at the anastomotic site.
- Prior to completion of the anastomosis, just before tying the final

sutures, gently remove the shunt.

WARNINGS

- Use care during insertion and removal of the shunt to avoid tearing the vessel and/or intimal dissection. NEVER force the shunt into the vessel.
- Carefully choose the shunt by matching the correct bulb diameter required for the internal diameter of the vessel.
- During the anastomosis of the graft to the artery, care must be exercised so the suture needle will not pierce the shunt. The silicone material is soft and similar to tissue in consistency. Thus, it is not easy to feel the needle entering the shunt.
- Once the shunt has been removed, do not reinsert as residual blood may remain in the shunt's lumen.

- Do not re-sterilize. The shunt is a single use device and is not intended for reuse.

DISCLAIMER OF WARRANTIES

Synovis Surgical Innovations (SSI), a division of Synovis Life Technologies, Inc., warrants that reasonable care has been used in the manufacture of this device. This warranty is exclusive and in lieu of all other warranties whether expressed, implied, written or oral, including, but not limited to, any implied warranties of merchantability or fitness. As a result of biological differences in individuals, no product is 100% effective under all circumstances. Because of this fact and since SSI has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after it leaves its possession, SSI does not warrant either a good effect or against an ill effect following its use. The manufacturer shall not be liable for any incidental or consequential loss, damage or expense arising directly or indirectly from the use of this device. SSI will replace any device which is defective at the time of shipment. No representative of SSI may change any of the foregoing or assume any additional liability or responsibility in connection with this device.

Rx ONLY CAUTION

Federal (USA) law restricts this device to sale by, or on the order of a physician.

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Flo-Thru Intraluminal Shunt[®] is a registered trademark of Synovis Life Technologies, Inc.
Flo-Thru Intraluminal Shunt[®] ist ein Warenzeichen von Synovis Life Technologies, Inc.
Flo-Thru Intraluminal Shunt[®] est un nom commercial appartenant à Synovis Life Technologies, Inc.
Flo-Thru Intraluminal Shunt[®] è un marchio di fabbrica della Synovis Life Technologies, Inc.
Flo-Thru Intraluminal Shunt[®] es una marca comercial de Synovis Life Technologies, Inc.
Flo-Thru Intraluminal Shunt[®] is een handelsmerk van Synovis Life Technologies, Inc.