
Vascular Probe

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SYMBOL DEFINITIONS:

	Ethylene Oxide sterilized
	Do not reuse
	Use by
	Lot
	Catalog number
	Consult <i>Instructions for Use</i> .
	Shaft length
	Bulb diameters

PN: SSI part number. TN: SSI tracking number. IC: SSI internal code.

INDICATIONS:

The disposable Vascular Probe is intended to be used during coronary and peripheral vascular surgery to probe blood vessels distally and proximally for blockages, to measure the internal diameter of vessels, and to act as an intravascular retractor.

DESCRIPTION:

The Vascular Probe is a sterile, single use, disposable polymeric device with a streamline bulb on each end of a flexible shaft. The bulbs are of different size on each end of the Vascular Probe. An appropriately-sized bulb is inserted either proximally or distally through the arteriotomy to probe the interior of the vessel and may be used by the surgeon to retract the vascular wall.

INSTRUCTIONS FOR USE:

Each surgeon is responsible for the proper procedure and technique to be used with this device.

The Vascular Probe is to be chosen based upon the internal diameter of the vessel. Never force the Vascular Probe into the vessel because sub-intimal dissection could occur. In addition to the exploration of the interior of the vessel and the detection of strictures, the device may also be used as a gentle intravascular retractor during completion of the anastomosis. Upon withdrawal from the distal and proximal portions of the anastomosis, the passage of the streamlined bulb confirms the patency and size of the just-created communication. A Vascular Probe of appropriate size may also be used as an occluder during surgery to allow maintenance of a dry surgical field.

WARNINGS:

- Because of some vessels' friability and the possible danger of tearing, care must be used when dealing with vessels.
- Guard against subintimal dissection. NEVER force the Vascular Probe's bulb into a vessel.
- Carefully choose the correct size the Vascular Probe according to the bulb diameter required.
- Do not resterilize. The Vascular Probe is a disposable, single-use device.
- Do not subject the Vascular Probe to the following agents: alcohol esters and ethers, any derivatives of phenols, compounds of mercury, chloride, peroxide, acetone, or any type of solvent.
- Do not use mechanical gripping or clamping devices to manipulate the probe.

DISCLAIMER OF WARRANTIES:


Synovis Surgical Innovations (SSI), a division of Synovis Life Technologies, Inc., warrants that reasonable care has been used in the manufacture of this device. This warranty is exclusive and in lieu of all other warranties whether expressed, implied, written or oral, including, but not limited to, any implied warranties of merchantability or fitness. As a result of biological differences in individuals, no product is 100% effective under all circumstances. Because of this fact and since SSI has no control over the

conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after it leaves its possession, SSI does not warrant either a good effect or against an ill effect following its use. The manufacturer shall not be liable for any incidental or consequential loss, damage or expense arising directly or indirectly from the use of this device. SSI will replace any device which is defective at the time of shipment. No representative of SSI may change any of the foregoing or assume any additional liability or responsibility in connection with this device.

Rx ONLY CAUTION:

Federal (USA) law restricts this device to sale by, or on the order of, a physician.

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AR-MED Limited (Company Number 3154742)
A Division of Synovis Life Technologies, Inc.